BD PureWick™ Flex
Female External Catheter
Instructions for Use

en	English	Female External Catheter	**2**
fr	French/Français	Sonde externe pour femmes	**6**
de	German/Deutsch	Externer Katheter für Frauen	**10**
it	Italian/Italiano	Catetere femminile esterno	**14**
es	Spanish/Español	Catéter externo para mujeres	**18**
nl	Dutch/Nederlands	Vrouwelijke externe katheter	**22**
pt	Portuguese/Português	Cateter Externo para Mulher	**26**
el	Greek/Ελληνικά	Εξωτερικός καθετήρας για γυναίκες	**30**
da	Danish/Dansk	Eksternt kateter til kvinder	**34**
sv	Swedish/Svenska	Extern kateter för kvinnor	**38**
fi	Finnish/Suomi	Naisten ulkoinen katetri	**42**
no	Norwegian/Norsk	Eksternt kateter for kvinner	**46**
pl	Polish/Polski	Cewnik zewnętrzny żeński	**50**
hu	Hungarian/Magyar	Női külső katéter	**54**
cs	Czech/Česky	Externí katétr pro ženy	**58**
ro	Romanian/Română	Cateter extern pentru femei	**62**
sk	Slovak/Slovenčina	Ženský externý katéter	**66**
tr	Turkish/Türkçe	Dişi Eksternal Kateter	**70**
ru	Russian/Русский	Женский наружный катетер	**74**
ko	Korean/한국어	여성용 외부 카테터	**78**
82	قسطرة خارجية أنثوية	العربية/Arabic	**ar**
86	کاتتر خارجی زنانه	فارسی (فارسی)/Persian (Farsi)	**fa**

BD PureWick™ Flex

en

Female External Catheter

Instructions for Use

Description:

The PureWick™ Flex Female External Catheter is made of three main parts, the white wicking material, the tubing, and the outer shell. The white wicking material absorbs voided urine and allows it to be collected at the bottom end of the device. The tubing draws the collected urine from the bottom end of the device to a collection canister and the purple shell holds the device together.

Indications

The PureWick™ Flex Female External Catheter is intended for non-invasive urine output management in adult users with female anatomy, for conditions such as urinary incontinence.

Intended users:

* User/Patient

* Caregiver/Healthcare professional


Contraindications

Users with urinary retention.


Warnings

- The PureWick™ Flex Female External Catheter is designed for single-use only. Trying to wash, sanitize or reuse it may lead to risk of infection or illness. This may damage the product and materials. The product may not work and cause injury or illness.
- To avoid potential skin injury, never push or rub the product against the skin during placement or removal.
- Do not block airflow to the product (e.g., blocking the vent hole, bedpan, barrier creams).
- Do not insert the product into vagina, anus, or other body orifice.
- Stop use if an allergic reaction occurs.
- After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state, and federal laws and regulations.
- Do not use on users with frequent episodes of stool incontinence without a fecal management system in place.
- Do not use on users with skin irritation or breakdown at the site.
- Use with caution on users who had recent surgery of the external female anatomy.
- Do not use on users with moderate/heavy menstruation who cannot use a tampon.

THE FEMALE ANATOMY

THE PRODUCT

MRI Warnings

- **Always disconnect the PureWick™ Flex Female External Catheter from wall suction or the PureWick™ Urine Collection System prior to an MRI procedure.**
- **It is important to closely follow these specific conditions that have been determined to permit the examination to be conducted safely. Any deviation may result in injury to the user.**
- Non-clinical testing demonstrated that the PureWick™ Flex Female External Catheter is MR Conditional. A user with one of these PureWick™ Flex Female External Catheters can be scanned safely under the following conditions.
- Static magnetic field of 1.5-Tesla or 3-Tesla only.
- Maximum spatial gradient magnetic field of 40-T/m (4000-Gauss/cm).
- Maximum MR system reported whole body averaged specific absorption rate (SAR) of 2-W/kg for up to 60 minutes of continuous scanning in the Normal Operating Mode.
- Under the scan conditions defined, the PureWick™ Flex Female External Catheter is expected to produce a maximum temperature rise of 1.5 degrees Celsius after 15-minutes of continuous scanning.
- There are no circularly polarized (quadrature-driven) RF coil restrictions.

Precautions

- Use with caution on users who are agitated, combative or uncooperative and might remove the product.
- Do not use barrier cream on the surfaces where the product will be placed. Barrier cream may reduce suction.
- Keep suction on while removing the product. This helps prevent leaks during removal.

Recommendations

- Check position of the product after repositioning the user.
- Users may transition (e.g., bed to chair), but the product should be removed if they are up and walking.
- Reposition the product and check the user's skin at least every two hours or per healthcare professional/facility recommendation.
- Placing the product snugly between the inner labia and buttocks holds it in place for most users. Breathable underwear may be useful for securing the product for some patients.
- Change the suction tubing and canister per facility protocol or per PureWick™ Urine Collection System Instructions for Use.
- Before connecting the product to wall suction tubing, verify suction function by covering the open end of the suction tubing with one hand and looking at the pressure dial. If the pressure does not increase when the line is covered, make sure that the tubing is secured, connected, and not kinked.

Instructions for Use

Wash hands before and after this procedure. If placing or removing this product for another person, wear gloves.

Setup

1 **If using continuous wall suction**, always use the minimum amount of suction necessary. Connect the canister setup per facility protocol. Set suction to 40 mmHg and adjust as needed. Different setups require higher suction typically between 60-100 mmHg.

NOTE: *If facility protocol allows and if using a graduated canister, captured urine may be used for approximate urine output measurement.*

If using the PureWick™ Urine Collection System, make sure all tubing connections are snug, and **turn it on**. Pressure adjustments are not necessary. Consult the PureWick™ Urine Collection System Instructions for Use as needed.

NOTE: *Ensure all connections are air-tight and check for possible cracks, leaks, kinks, or occlusions.*

2 Before placing the product, curve it to fit the user's anatomy. This helps the product stay in place.

3 **If using continuous wall suction**, securely connect the product to the suction tubing.

NOTE: *Keep track of how long the product has been in place by writing down the date and time it was placed.*

If using the PureWick™ Urine Collection System, securely connect the PureWick™ Flex Female External Catheter to the collector tubing.

Placement

4 Have the user lie on their back with knees bent and thighs apart (frog legged) or on their side before placing the product.

With their legs open, check the skin for irritation or breakdown, and clean the female anatomy.

NOTE: *If skin irritation or breakdown is seen, do not use the product.*

5 Spread the inner labia and keep them spread while placing the product.

With white wicking material side facing the user, place the bottom end of the product between the buttocks, but not as far back as to cover the anus.


Gently tuck white wicking material side between the inner labia, directly against the urethra (where the urine comes out).


NOTE: *Make sure the urethra is at least one inch down from the top of the white wicking material.*


6 Make sure the inner labia wrap **around** the product and are not tucked under it.


Once the product is in place, slowly close the legs and **make sure the vent hole is not covered**. See product diagram on page 1 for vent hole location.


NOTE: *For users with **larger** labia, less of the product will be visible. For users with **smaller** labia, more of the product will be visible.*



Removal

7 Open the legs and gently spread the labia. Gently pull the product away from the body. Keep suction on while removing the product.

After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state, and federal laws and regulations.

Maintenance:

8 Replace the product every 12 hours or if soiled with feces or blood.

SYMBOL AND SYMBOL TITLE (STANDARD REFERENCE)			
Symbols immediately below are derived from the following standard: ISO 15223-1:2021 - Medical devices — Symbols to be used with medical device labels, labeling, and information to be supplied – Part 1: General requirements.			
	Manufacturer	eifu.bd.com	Consult instructions for use or consult electronic instructions for use
LOT	Lot number		Importer
REF	Catalog number	EU REP	Authorized representative in the European Community
	Use by	CH REP	Authorized representative in Switzerland
NON STERILE	Non-sterile	UK REP	Authorized representative in the United Kingdom
	Do not use if package is damaged and consult instructions for use	MD	Medical Device
	Single use	UDI	Unique Device Identifier
Symbols immediately below are derived from the following standard: ISO 7000 – Graphical symbols for use on equipment			
	Packaging unit		
Symbols immediately below are derived from the following standard: ASTM F2503 – Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment			
MR	Magnetic Resonance (MR) Conditional		
Other			
	Not made with Natural Rubber Latex		

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BD, the BD Logo and PureWick are trademarks of Becton, Dickinson and Company or its affiliates. © 2025 BD. All rights reserved.
C. R. Bard, Inc.
EU REP Becton Dickinson Ireland Ltd.
Becton Dickinson AG, Biningerstrasse 94
UK REP Becton, Dickinson U.K. LTD